

Lens Parameters

Bausch + Lomb ULTRA® ONE DAY Multifocal

MATERIAL:	kalifilcon A
LENS MATERIAL TECHNOLOGY:	Advanced MoistureSeal® Technology + ComfortFeel Technology
WATER CONTENT:	55%
OXYGEN TRANSMISSION:	134 Dk/t @ -3.00D
LENS DESIGN TECHNOLOGY:	3-Zone Progressive™ Design
BASE CURVE:	8.6 mm
DIAMETER:	14.2 mm
CENTRE THICKNESS:	0.08 mm @ -3.00D
SPHERICAL POWERS:	+6.00D to -10.00D in 0.25D
ADD POWERS:	Low: up to +1.50D, High: up to +2.50D
VISIBILITY TINT:	Light blue
MODALITY:	Daily
UV PROTECTION:	Yes (CLASS II)
MODULUS:	0.5MPa
PACK SIZE:	30 and 90 pack

References: † when the fitting guide was followed 1. Results from a 20-site, 3-week study of Bausch + Lomb (kalifilcon A) Daily Disposable Multifocal contact lenses on 294 habitual multifocal soft contact lens wearers. Only Bausch + Lomb ULTRA® ONE DAY Multifocal contact lenses combine a 3-Zone Progressive™ Design for effortless clarity from near to far and a complete moisture + comfort system with Advanced MoistureSeal® and ComfortFeel Technologies that offer a moisture rich environment, high Dk/t, low modulus, UV blocking. Bausch + Lomb ULTRA® ONE DAY Multifocal contact lenses deliver health through its complete system working together to support a healthy ocular surface environment, the inclusion of eye health ingredients which are retained over 16 hours and the high allowance of oxygen permeability (Dk/t=134).

Bausch + Lomb ULTRA® ONE DAY, MoistureSeal, and 3-Zone Progressive are trademarks of Bausch + Lomb Incorporated or its affiliates. Bausch + Lomb Canada.

®/™ are trademarks of Bausch + Lomb Incorporated or its affiliates. Bausch + Lomb Corporation, Vaughan, Ontario, L4K 4B4 © 2024 Bausch + Lomb. TP21150 160-8892E

BAUSCH + LOMB
See better. Live better.

FITTING GUIDE FOR Bausch + Lomb ULTRA® ONE DAY Multifocal

STEP 1: Update spectacle refraction and Add power

STEP 2: Select contact lens distance prescription based upon spherical equivalent from spectacle Rx and following Add guidance (adjusted for vertex distance if necessary)

ADD SELECTION:

SPECTACLE Add	BOTH EYES
+0.75D to +1.50D	Low Add
+1.75D to +2.50D	High Add

EVALUATE THE LENS FOR SUCCESS

- Allow trial lenses to equilibrate for at least 10 minutes before assessing fit and vision
- Evaluate distance and near vision binocularly in normal room illumination
- If vision at distance and near are satisfactory, dispense lenses and schedule follow-up exam within 1-2 weeks



REFINE IF NEEDED

Determine eye dominance at distance by placing a +1.50 loose hand-held trial lens alternately over each eye binocularly through updated distance correction. The eye for which binocular vision is blurriest through the +1.50 is the dominant eye.

Peers demonstrated success with the 3-Zone Progressive™ Design†1:

- Easy to fit for 99% of patients
- 87% of patients successfully fit in one visit
- 99% of patients successfully fit in two visits

		Near Vision		Distance Vision				
If patient is wearing:	TWO LOW ADDS	DOMINANT EYE		NON-DOMINANT EYE				
		Initial Lens	Low Add	Low Add		Initial Lens	Low Add	Low Add
		Refinement 1	Low Add	High Add		Refinement 1	Bausch + Lomb ULTRA® ONE DAY Sphere	Low Add
		Refinement 2: If vision is still unsatisfactory, make small changes by adding +0.25D at a time to non-dominant eye (wearing High Add lens) using hand-held lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.				Refinement 2: If vision is still unsatisfactory, make small changes by adding -0.25D at a time to dominant eye (wearing Bausch + Lomb ULTRA® ONE DAY spherical lens) using hand-held lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.		
	TWO HIGH ADDS	DOMINANT EYE		NON-DOMINANT EYE				
		Initial Lens	High Add	High Add		Initial Lens	High Add	High Add
		Refinement 1	High Add	Add +0.25D to the non-dominant eye		Refinement 1	Low Add	High Add
		Refinement 2: If vision is still unsatisfactory, make small changes by adding +0.25D at a time to non-dominant eye using hand-held lenses, and continue evaluating vision binocularly at normal room illumination. Adjust contact lens power when vision is satisfactory.				Refinement 2: If vision is still unsatisfactory, make small changes by adding -0.25D at a time to dominant eye (wearing Low Add lens) using hand-held lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.		

Scan to use the Bausch + Lomb
Multifocal Fitting calculator!

